

During Zepto capsulotomy, most OVD with a viscosity of $\leq 300,000$ mPas may be used. Examples:	Viscosity (mPas)
Cohesive Healon Pro Provisc Amvisc Amvisc Plus NuVisc	150 K 280 K 100 K 60 K 45 K
Viscous dispersive DisCoVisc	100 K
Dispersive Healon Endocoat Viscoat Cellugel Vitrax	50 K 50 K 40 K 25 K
Dispersive very low viscosity i-Cell Ocuvis Hymecel Adatocel Visilon	6 K 4.3 K 4 K 4 K 4 K

After Zepto capsulotomy, any OVD may be used. Examples:	Viscosity (mPas)
Viscoadaptives Healon 5 Pro Healon GV Pro MicroVisc (iVisc) Plus	4 M 2 M 3.3 M
Viscous-cohesive viscoelastic MicroVisc (iVisc) Phaco MicroVisc (iVisc)	1 M 1 M
Dispersive ClearVisc OcuCoat	40 K 4 K